

OptraSculpt®

Efficient contouring of composites

Time saving
contouring due to
non-stick effect



OptraSculpt®

Efficient contouring of composites

The contouring and shaping of state-of-the-art high-quality dental composites is a challenging task. Apart from other factors, it significantly influences the long-term success of the treatment, since the correct anatomical shape of the restoration is crucial to the overall health of the chewing system.

The occlusal design is important for the proper functioning of the masticatory system. Correctly designed proximal marginal ridges facilitate the cleaning of interdental spaces. Moreover, smooth restoration surfaces minimize the accumulation of plaque, and carefully created margins reduce the occurrence of secondary caries.

Innovative solutions for the dental practice

**Also suitable for contouring
lab composites**

Easy adaptation

- **Non-stick properties:** reduce the stickiness of the composite versus the instrument
- **Compatibility:** attachments are suitable for both instruments
- **Working angles are adjusted** to the anterior and posterior regions

High efficiency

- Efficient, easy contouring due to specially adjusted attachment shapes
- Reduced effort in finishing and polishing

Professional esthetics

- No marks left by the instrument on the restoration surface
- Reduction of air inclusions in the composite during application
- Homogeneous surface texture



OptraSculpt® Pad

for anterior teeth

- **Foam pad attachments** for non-stick, mark-free contouring
- **Two pad sizes** for efficient contouring of direct veneers and large-surface class III and IV restorations
- **Reference scale** on the instrument handle for recording the axis alignment and the width of the anterior teeth
- **Optimum working angle** for the anterior region




OptraSculpt® Next Generation

for posterior teeth

NEW

- **Three attachment shapes** for the professional creation of sophisticated anatomical tooth structures: **ball, pointed tip, chisel**
- **Optimum working angle** for the posterior region
- **Hygienic storage container** for easy access to attachments




Clinical application

OptraSculpt® Pad

Contouring in comparison



Easy and mark-free with OptraSculpt® Pad



Complicated and time-consuming with a metal spatula

Orientation aids on the instrument handle




Reference scale 1 for estimating the tooth width



Reference scale 2 for determining the tooth axis

OptraSculpt® Next Generation



Ball – adaptation of composites

- Box floor
- Deep increments



Pointed Tip – professional contouring of fissures

- Fissures
- Cusp slopes



Chisel – creation of proximal marginal ridges

- Proximal marginal ridges
- Cusp slopes
- Class V cavities

Delivery forms

OptraSculpt® Next Generation & OptraSculpt® Pad

683066AN	OptraSculpt & OptraSculpt Pad System Kit 1 x single-ended instrument 1 x double-ended instrument 60x pads / 4 mm 20x ball 20x chisel 20x pointed tip
----------	---

OptraSculpt® Pad

639791AN	OptraSculpt Pad Assortment 1 x single-ended instrument 60x pads / 4 mm 60x pads / 6 mm
639792AN	OptraSculpt Pad Refill / 4 mm / 60 60x pads / 4 mm
639793AN	OptraSculpt Pad Refill / 6 mm / 60 60x pads / 6 mm
646004AN	OptraSculpt Pad Instrument Refill 1 x single-ended instrument

NEW

OptraSculpt® Next Generation

683067AN	OptraSculpt Starter Kit 1 x double-ended instrument 20x ball 20x chisel 20x pointed tip
683068AN	OptraSculpt Refill / ball / 60 60x ball
683069AN	OptraSculpt Refill / chisel / 60 60x chisel
683100AN	OptraSculpt Refill / point / 60 60x pointed tip
683154AN	OptraSculpt Instrument Refill 1 x double-ended instrument



Direct Restoratives

OptraSculpt® forms a part of the "Direct Restoratives" product category. The products of this category cover the procedure involved in the direct restoration of teeth – from preparation to restoration care. The products are optimally coordinated with each other and enable successful processing and application.



THESE ARE FURTHER PRODUCTS OF THIS CATEGORY:

Tetric EvoCeram® Bulk Fill & Tetric EvoFlow® Bulk Fill

High-performance posterior composite



The efficient posterior composite

- Increment thickness of up to 4 mm due to the highly reactive light initiator Ivocerin®
- Dentin- or enamel-like volume replacement
- 10 seconds ($\geq 1,000$ mW/cm²)
- 47% less time required than with the conventional technique*

* compared to Tetric EvoFlow® and Tetric EvoCeram®. Data available on request.

Bluephase® Style

The curing light



The smallest LED for every use

- Every hand – comfortable to hold for men and women
- Every material – universal use due to Polywave® LED with broadband spectrum
- Every time – optional corded operation due to Click & Cure

Would you like to know more about the products of the "Direct Restoratives" category? Simply get in touch with your contact person at Ivoclar Vivadent or visit www.ivoclarvivadent.com for more information.

Ivoclar Vivadent AG
Bendererstr. 2
9494 Schaan
Liechtenstein
Tel. +423 235 35 35
Fax +423 235 33 60
www.ivoclarvivadent.com

Descriptions and data constitute
no warranty of attributes.

Printed in Germany
© Ivoclar Vivadent AG, Schaan/Liechtenstein
684580/en/2016-09

ivoclar
vivadent
passion vision innovation