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Chrom-Kobalt-Bonding ckb

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Instructions for use

Chrome-Cobalt-Bonding ckb

Before using this product, please read these instructions for carefully use.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Description of the product

Chrome-cobalt bonding (ckb) is a ceramic material. It is used to form a layer between all metals and ceramic. It compensates for different expansion coefficients between metal and ceramic and blocks escaping metal oxides.

2. Indication

Chrome-cobalt bonding (ckb) is used for non-precious metal alloys, titanium alloys, precious metal and reduced precious metal content alloys.

3. Contraindication

So far the following alloys have proved to be unsuitable: Crutanium by Krupp Austenal, Crystalloy by Shofu and Ticonium.

4. Safety recommendations and hazard warnings

4.1 Protective clothing

While processing chrome-cobalt bonding, dust mask, safety gloves, goggles and protective clothing must be worn.

4.2 To ensure safe handling of the product, the formation of dust must be avoided.

4.3 After skin contact: Immediately wash off with water and soap and rinse thoroughly. Normally, the product does not irritate the skin.

After eye contact: Rinse eyes carefully with running water with the lid being opened.

After ingesting: If problems persist, seek medical attention.

5. Storage and durability

5.1 The product must be sealed well and stored dry. If proper storage is ensured, the product has almost unlimited shelf life.

5.2 Proper disposal

Product can be disposed of with household garbage.

6. Processing

6.1 Preparation of the object

After devesting clean the surface of the framework with a pencil blaster using aluminium oxide 110µm. Clean it carefully with distilled water in the ultrasonic unit or with a steam-blasting unit. After this process, the metal surface must not be touched any more.

6.2 Mixing

Mix the chrome-cobalt bonding with distilled water to obtain a thin consistency. Other mixing liquids are not suitable.

6.3 Applying the chrome-cobalt bonding

Use a staining brush to apply a single covering layer of the mixed chrome-cobalt bonding onto the metal surface. A uniform, pink color should result on the surface.

Please note: Once chrome-cobalt bonding has dried, it must not be mixed anew!

6.4 Firing of the chrome-cobalt bonding in the ceramic furnace

Firing of chrome-cobalt bonding is started at 650° C and is continued up to a temperature of 980° C. The rate is 55° C/minute. The firing time amounts to 6 minutes. Firing is performed under vacuum. After firing, the framework must be immediately taken out of the furnace. The framework should have a golden-yellow to beige color after firing.

6.5 Opaque firing and further processing

A covering layer of opaque - which corresponds to the tooth shade - is applied onto the fired chrome-cobalt bonding. A washbake is not required. Then continue according to the instructions of the manufacturer of the ceramic material.

6.6 Sources of errors and how to eliminate them

Since the middle of 1982 chrome-cobalt bonding has been used every day in numerous laboratories. To ensure successful application we recommend to test the desired alloy with chrome-cobalt bonding.

6.6.1 Cracks in the chrome-cobalt bonding: Too much chrome-cobalt bonding has been applied! Remove the material by sandblasting and apply anew as described under item 6.3.

6.6.2 Bubbles in the chrome-cobalt bonding: Degassing of the metal is required. Remove chrome-cobalt bonding by sandblasting and - prior to firing anew - degas the framework under vacuum at 980° C.

6.6.3 Green and/or black stains in the chrome-cobalt bonding: chrome-cobalt bonding was mixed too thinly or the framework was moist or greasy. Stains may also result if dried chrome-cobalt bonding has been mixed anew. Remove chrome-cobalt bonding by sandblasting, clean the framework and then apply chrome-cobalt bonding anew.

6.6.4 Cracks in the ceramic: If chrome-cobalt bonding displays a rust-like brown color after firing, the alloy is not suitable. Please use a different alloy. Up until today the following alloys have proved to be not suitable: Crutanium by Krupp Austenal, Crystalloy by Shofu and Ticonium. Cracks may also result, if dried chrome-cobalt bonding has been mixed anew since it is not ensured that the different expansion coefficients are compensated.

7. Additional information

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on processing. Any occurring claims for damages may only be made up to the value of our products.

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