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Safety instructions and general instructions for BEGO investment materials

en

Safety instructions

- Investment materials contain quartz. **Do not inhale dust!** Risk of damage to lungs (silicosis, lung cancer).
Recommendation: Use gas mask type FFP 2 – EN 149:2001. Open bag along the perforation or cut open with scissors and avoid generation of dust when filling into mixing bowl. Rinse out empty bag with water before crumpling up. **Following eye contact:** Remove contact lenses, hold eyelids open and rinse thoroughly with clean, running water for at least 15 minutes. Consult an eye specialist.
 - Remove dust at workplace only when **moist**.
 - To avoid dust during deflasking, place the mould in water after it has completely cooled down after casting until it is thoroughly moistened.
 - Use suction extraction system with fine dust filter when blasting.
- Remove dried-up mixing liquid only when **moist** (fine dust). Seal bottle securely after use.
- The setting reaction of the investment-material generates heat. To avoid injuries, please use gloves or suitable muffle tongs.
- Combustion gases resulting from the preheating of casting moulds must be discharged via a suitable smoke outlet.
- Only with shock heat investment materials: **Risk of injury during shock heating!** Place all moulds in the furnace within 10 seconds – then keep the furnace door closed for 15 minutes!
- Never touch hot moulds with your hands. Use suitable mould tongs to avoid burns.
- After casting, allow the moulds to cool down until warm to the touch in a protected location!

Danger

Contains: Quartz (SiO₂), Cristobalite
H372i Causes damage to organs through prolonged or repeated exposure if inhaled.
P260 Do not breathe dust.
P314 Get medical advice/attention if you feel unwell.
P501 Dispose of contents/container to a facility in accordance with local and national regulations.

General instructions

- The ideal processing temperature for BEGO investment materials is **20 °C** (except Bellasun: 30 °C). To keep it constant at a higher ambient temperature, use a temperature control cabinet if necessary. At a high ambient temperature also place the mixing bowl and mixing bowl in the temperature control cabinet.
- Observe shelf life dates of investment material and mixing liquid! Do not use without checking after the specified date. Do not bring phosphate-bonded investment materials into contact with plaster or investment material containing plaster.
- Do not use crystallized mixing liquid any longer. The higher the concentration of mixing liquid, the greater the expansion. **A practical tip: Keep spare bottles!** Fill BegoSol® bottle up to desired % mark with BegoSol® and top up to 100 % with distilled water. Label this spare bottle with % data.
Applies accordingly to BegoSol HE/K.
Advantage: The mixing liquid is available with the desired concentration and can be measured in one operation.
- Data on the concentration of the mixing liquid are based on our own experience and tests and are regarded as reference values.
- Packaging must be completely emptied and disposed of properly in accordance with the statutory regulations. Rinse out the empty bag with water before compressing.

Application-related recommendations made by us, whether given verbally, in writing or by way of practical instructions, are based on our own experience and tests and may therefore be regarded only as general guidelines.

Our products are subject to continuous development. We therefore reserve the right to change the design and composition of our products. Work results may be influenced by such parameters as temperature, mixing intensity, investment system, etc.