ORBIS Cavity Liner LC

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PRODUCT DESCRIPTION:

Orbis Cavity Liner LC is a light-curing, radiopaque one-component material containing calcium hydroxide.

FIELDS OF APPLICATION:

Indirect pulp capping
As a cavity liner under all restorative materials

APPLICATION:

During application, ensure that the tooth surfaces are perfectly dry and avoid contamination. Instead of using **Orbis Cavity Liner LC** use, an aqueous calcium hydroxide product on dentine close to the pulp (< 0.5 mm) and for direct pulp capping.

Orbis Cavity Liner LC can be applied directly from the syringe. Dispense the required amount of **Orbis Cavity Liner LC** from the tube onto a mixing pad and apply using a suitable instrument.

Cavity linings with a thickness in excess of 1 mm need to be built up in layers and each layer cured using blue light for at least 20 s.

If increased adhesion values are required (e. g. when using **Orbis Cavity Liner LC** for adhesive restorative treatments), a suitable adhesive can be used prior to the application of **Orbis Cavity Liner LC**. To this end, follow the adhesive manufacturer's instructions for use.

Carefully remove any excess material prior to light-curing using suitable tools. Do not remove polymerised excess material with a probe but, instead, carefully use a rotating instrument and then place the restoration according to the manufacturer's instructions.

INFORMATION, PRECAUTIONARY MEASURES:

Do not use **Orbis Cavity Liner LC** in the restoration margins. In the marginal area, the permanent restoration should always be placed on caries-free dental hard tissue which, when using the adhesive technique, must have been conditioned accordingly.

Orbis Cavity Liner LC adheres to methacrylate-based products, e. g. crown and bridge materials, temporary sealing materials and certain impression materials (e. g. polyether), which can lead to removal of the cavity lining. In such cases, application of a separating agent is recommended (e. g. Vaseline, glycerine gel etc.). Remove the insulating material prior to placement of the final restoration.

Once dispensed, **Orbis Cavity Liner LC** must not be exposed to direct light (operating light, sunlight) and must be applied as quickly as possible.

Orbis Cavity Liner LC contains methacrylates, BHT and amines. **Orbis Cavity Liner LC** should not be used in case of known hypersensitivities (allergies) to any of these ingredients. Calcium hydroxide has a caustic effect on the skin and mucosal tissue, it must not be allowed to come into contact with the eyes!

Phenolic substances, especially products containing eugenol and thymol, impede the curing of **Orbis Cavity Liner LC**. The use of zinc oxide eugenol cements or other eugenol-containing materials in combination with **Orbis Cavity Liner LC** should therefore be avoided.

Hold the output window of the light as close to the surface of the restoration as possible, as curing may otherwise be insufficient. Inadequate curing may lead to loss of strength and to discomfort.

STORAGE:

Store at temperatures of 4° C to 23° C protected from light and moisture.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and / or written information has been given to the best of our knowledge and without obligation. Our information and / or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.