DENTOCORE Body

CORE BUILD-UP

DENTOCORE Body

CORE BUILD-UP

DENTOCORE Body

CORE BUILD-UP

DentoCore Body Core Build-up

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DentoCore Body Core Build-up

Nanofill composite core build-up material. Dual-cure.

INDICATIONS:

DentoCore Body is a dual-cure, fluoride releasing, resin composite material designed for the fabrication of core build-ups. DentoCore Body may be used with any Bis-GMA compatible bonding agent. DENTOCORE and DENTOCORE Body having the same properties and indications, the choice of the material depends on dentist habits.

PROPERTIES:

• Easy to use. • Provides a consistent, homogeneous mix. • Dual curing: material can be light cured, additional chemical cure assures that the material can be quickly, bulk filled. • Dense void-free core cuts like dentin. • Fluoride release: strengthens natural tooth structure and helps prevent secondary caries. • Radiopaque. • Superior compressive strength for long-lasting restorations. • Enhanced mechanical properties of conventional Bis-GMA composite by utilization of Interlaced and Exfoliated nanocomposite technology.

INSTRUCTIONS FOR USE:

1) Prepare the tooth according to clinical requirements. 2) Isolate the tooth from oral fluids. A rubber dam is recommended. 3) Acid etch the entire preparation with a 32%-40% phosphoric acid for 15 seconds or apply a self etching adhesive (e.g. Quick Bond of Itena) and follow instructions for use of this manufacturer. 4) Rinse off etchant with a water spray for at least 15 seconds and apply bonding agent of choice (e.g. DentoBond Adhesive of Itena). 5) Place into base of prepared tooth an adequate amount of DentoCore Body to complete build-up. 6) Do not disturb DentoCore Body during the gelation stage. Light-cure for at least 40 seconds. Check for complete hardness. Chemical-cure intraoral set-time is approximately 4 1/2 minutes from start of mix/application. Settime will be longer outside of mouth.

TECHNICAL DATA:

Working Time:	1:30 minute
Setting Time:	4:30 minutes
Compressive strength:	250 MPa
Diametral tensile strength:	40 MPa
Flexural strength:	200 MPa
E-Modulus:	8 GPa
Water sorption:	14 µg/mm³
Solubility:	2 µg/mm³

CAUTION!

DentoCore Body contains polymerisable monomers which may cause allergic contact dermatitis in susceptible persons. Wash thoroughly with soap and water after contact. If skin sensitization occurs, or if known allergy to methacrylate resin exists, discontinue use.

PRECAUTIONS IN USE:

• Do not use in presence of eugenol based materials; complete setting of DentoCore Body may not take place. • A calcium hydroxide or glass ionomer liner may be used to isolate DentoCore Body from eugenol sources. • If refrigerated, allow to return to room temperature before use. Refrigeration will increase viscosity, making application more difficult, and will retard setting time.

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CORE BUILD-UP

DentoCore Core Build-up

EN
Flowable nanofill composite core build-up material. Dual-cure.

INDICATIONS:

DentoCore is a dual-cure, fluoride releasing, flowable resin composite material designed for the fabrication of core build-ups. DentoCore may be used with any Bis-GMA compatible bonding agent.

PROPERTIES:

• Easy to use. • Provides a consistent, homogeneous mix. • Dual curing: material can be light cured, additional chemical cure assures that the material can be quickly, bulk filled. • Dense void-free core: cuts like dentin. • Fluoride release: strengthens natural tooth structure and helps prevent secondary caries. • Radiopaque. • Superior compressive strength for long-lasting restorations. • Enhanced mechanical properties of conventional Bis-GMA composite by utilization of Interlaced and Exfoliated nanoparticle technology.

INSTRUCTIONS FOR USE:

1) Prepare the tooth according to clinical requirements. 2) Isolate the tooth from oral fluids. A rubber dam is recommended. 3) Acid etch the entire preparation with a 32%-40% phosphoric acid for 15 seconds or apply a self etching adhesive (e.g. Quick Bond of Itena) and follow instructions for use of the manufacturer. 4) Rinse off etchant with a water spray for at least 15 seconds and apply bonding agent of choice (e.g. DentoBond Adhesive of Itena). 5) Place into base of prepared tooth on adequate amount of DentoCore to complete build-up. 6) Do not disturb DentoCore during the gelation stage. Light-cure for at least 40 seconds. Check for complete hardness. Chemical-cure intraoral set-time is approximately 4 1/2 minutes from start of mix/application. Set time will be longer outside of mouth.

TECHNICAL DATA:

Working Time:	1 : 30 minute
Setting Time:	4 : 30 minutes
Compressive strength:	250 MPa
Diametral tensile strength:	40 MPa
Flexural strength:	200 MPa
E-Modulus:	8 GPa
Water sorption:	14 µg/mm ³
Solubility:	2 µg/mm ³

CAUTION!

DentoCore contains polymerisable monomers which may cause allergic contact dermatitis in susceptible persons. Wash thoroughly with soap and water after contact. If skin sensitization occurs, or if known allergy to methacrylate resin exists, discontinue use.

PRECAUTIONS IN USE:

• Do not use in presence of eugenol based materials; complete setting of DentoCore may not take place. • A calcium hydroxide or glass ionomer liner may be used to isolate DentoCore from eugenol sources. • If refrigerated, allow to return to room temperature before use. Refrigeration will increase viscosity, making application more difficult, and will retard setting time.

ESSENTIAL INGREDIENTS:

• Bis-GMA • 2-Hydroxyethylmethacrylate • Triethyleneglycoldimethacrylate • Benzoyl peroxide • Co-initiators • Photoinitiators • Dendritic polymer • Barium aluminoborosilicate glass • Fluoride releasing filler • Fumed silica • Nano filler • White pigment, pink pigment or Shade A3-Alkyl quaternary ammonium bentonite.

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