

MT plus

Wet and dry trimmer

The MT plus is a wet and dry trimmer. Pure power, even the hardest stone is easily removed thanks to the powerful motor performance.

Advantages

- The front door can be easily removed without the use of tools to allow easy, quick maintenance.
- Detachable trimmer table and spray tube enable very easy cleaning.
- Noise-reducing, scratch-resistant housing made from a special material (polyurethane).

Details

- The trimmer can be used for any model situation due to infinitely adjustable trimmer table.
- Integrated tool allows quick change of the trimmer disc.
- The motor and water is stopped by a safety switch when opening the door, providing a high degree of safety.



tool in-
cluded



easy use



Powerful wet and dry trimming with minimum vibrations

Technical data

| | |
|---|---|
| Permissible mains voltage | 220–230 V |
| Permissible mains frequency | 50 Hz |
| Power consumption | 1300 W // 1.74 hp |
| Motor performance (P1) | 1300 W // 1.74 hp |
| Motor performance (P2) | 900 W // 1.21 hp |
| Rotational speed | 2850 rpm |
| Water consumption max. | 9 l/min // 0.32 cfm |
| Water pressure min. / max. | 1–5 bar // 14.5–72.5 psi |
| Ø (Connection piece for drain hose / external extraction) | 36 mm // 1.42 " |
| Ø (Trimmer disc) | 234 mm // 9.2 " |
| Dimensions (W x H x D) | 300 x 320 x 410 mm // 11.8 x 12.6 x 16.2" |
| Weight (without trimmer disc) | ~13,5 kg // ~30 lbs |

Ordering information

| | |
|--|--------------|
| MT plus incl. Klettfix trimmer disc, 220–230 V | No. 18030000 |
| MT plus incl. Marathon trimmer disc, 220–230 V | No. 18030500 |
| Trimmer discs | Page 63 |