



CanalPro™ CL (Cordless Handpiece)



- Ergonomic design and control
- Big, user friendly LCD area
- 5 storage programs
- Auto-reverse mod
- 6 head positions allow contra angle head to be adjusted to your preferred position
- Operating time: approx. 80 min.
- Quick charging time of approx. 90 min.
- Contactless charging via induction
- Power saving mode after 10 min. when not in use
- Contra angle head is autoclavable at 135 °C
- Rotation Speed: 140 ~ 550 rpm
- Torque Level: 0.3 ~ 3.0 Ncm

CanalPro CL & Accessories

Complete Set EU / UK (without Apex Locator connection) 1 × Motor Handpiece, 1 × Charger, 1 × Contra Angle 16:1	60011035 / 60011036	
Complete Set EU / UK (with Apex Locator connection) 1 × Motor Handpiece, 1 × Charger, 1 × Contra Angle 16:1 AL	60014050	
Motor Handpiece	1 pc	60011037
Charger	1 pc	60011038
Contra Angle 16:1 (without Apex Locator connection)	1 pc	60011039
Contra Angle 16:1 AL (with Apex Locator connection)	1 pc	60013747

CanalPro™ CL 2 (Cordless Handpiece)



- Ergonomic design and control
- Big, user friendly LCD area
- Multiple speed settings, Torque Level, Auto-reverse mod and Translation ratio
- Audio feedback for Auto-reverse mod
- Acoustic warning signal by 75 % of the set torque value
- 6 storage programs
- Operating time: approx. 90 min.
- Rotation Speed: 120 ~ 500 rpm
- Torque Level: 0.3 ~ 3.0 Ncm

CanalPro CL 2 & Accessories

Complete Set EU 1 × Motor Handpiece, 1 × Motor Handpiece Cover, 1 × Charger, 1 × Charging Cable, 1 × Contra Angle 16:1, 1 × Oil supply	60019362	
Motor Handpiece	1 pc	60018730
Motor Handpiece Cover	1 pc	60018735
Charger	1 pc	60018733
Charging Cable (EU)	1 pc	60019371
Contra Angle 16:1	1 pc	60018731
Oil supply	1 pc	60018736