

Tip changing made easy.

World premiere!

A single revolution

100 KaVo



Innovation and Tradition. Advantages for you:

- A faster and more comfortable tip change enables a smoother and more efficient procedure
- The same outstanding product features as the classic SONICflex 2003L
- Also works with Sirona couplings

New specialised tips for surgical indications

- Membrane Disc: Ideal for reliably and gently repositioning the Schneider membrane during sinus augmentation
- Saw: Ideal to allow precise, extremely narrow incisions, when cutting ligaments
- Bone Harvester: Ideal for easily and quickly harvesting large quantities of bone during sinus augmentation



KaVo. Dental Excellence.