



LC-TRAY PRO

 **dentalfarm**
DENTAL EQUIPMENT SINCE 1972

SIMPLY INCOMPARABLE

LC-TRAY PRO

SIMPLY INCOMPARABLE

FORNO DI POLIMERIZZATORE A LUCE LED PER PLACCHE

- FORNETTO SPECIFICO PER IL TRATTAMENTO DELLE PLACCHE IN RESINA FOTOSENSIBILE, PER LA PREPARAZIONE DI PORTA IMPRONTE INDIVIDUALI, BASI PER DENTIERE, MASTICONI DI PROVA, BLOCCHI OCCLUSALI E MASCHERINE
 - 36 HIGH-POWER UV LED - POTENZA 48W
 - BASE INTERNA RIFLETTENTE, CON AMPIO SPAZIO PER IL POSIZIONAMENTO DEGLI OGGETTI
 - PRATICO CASSETTO ESTRAIBILE FORNITO DI SCHERMO PROTETTIVO
 - DISPLAY LCD
 - TIMER ELETTRONICO. 4 PULSANTI PER LA SELEZIONE DI CICLI DA 3, 4, 5 E 6 MINUTI
-

LED LIGHT-CURING BOX FOR RESIN PLATES

- DESIGNED FOR THE PREPARATION OF INDIVIDUAL CUSTOM TRAYS, DENTURE BASEPLATES, TESTING BITE PLATES, FIXING MASKS AND OCCLUSAL BLOCKS
- FEATURES 36 HIGH-POWER UV LED - 48W
- UNIQUE DESIGN AND ENHANCED INTERNAL CAPACITY
- PRACTICAL SLIDING DRAWER WHERE PARTS CAN BE PLACED IN, FITTED WITH BLACK PROTECTION SHIELD
- LCD DISPLAY SCREEN
- ELECTRONIC TIMER PROVIDING FOUR DIFFERENT AVAILABLE TIME CYCLES OF 3, 4, 5 AND 6 MINUTES



A5413 LC-TRAY PRO

DIMENSIONI mm: 220 x 240 x 110h	<i>DIMENSIONS mm: 220 x 240 x 110h</i>
CASSETTO mm: 155 x 145 x 70h	<i>DRAWER mm: 155 x 145 x 70h</i>
ALIMENTAZIONE: 100/230V - 50/60Hz	<i>VOLTAGE: 100/230V - 50/60Hz</i>
ASSORBIMENTO: 1,5amp	<i>RATED POWER: 1,5amp</i>
36 HIGH-POWER LED - 48W	<i>36 HIGH- POWER LED - 48W</i>
SPETTRO: 365-405nm	<i>WAVELENGHT RANGE: 365-405nm</i>