

A6020 INDUCTION BURNER E-BUNSEN



Product group

Macchine elettriche per la lavorazione delle cere

Main features

- ⇒ Best safety thanks to no-flame technology
- ⇒ Compact size and possibility to position the equipment everywhere in the lab, provided an electrical plug is available
- ⇒ Positioning sensor to activate heating
- ⇒ High thermal yield, the tool will be hot in 2 seconds

Line tension

230v A.C. 50 Hz

Absorption

100 W

Length

90 mm

Width

210 mm

Height

90 mm

Net weight

0,9 kg

Gross weight

1,3 kg