

For a clean lab and dental practice



Grinding box

Acrylic grinding box with built-in magnifying glass prevents the spread of grinding dust.

- resistant, light box for flexible use within the lab.
- optimum control via a 2 x magnification magnifying glass.
- large volume and wide hand inserts.
- protective glasses and mouth masks are no longer required.

- ideal for working on small work pieces as nothing can fall to the ground or be sucked up by the extractor unit.
- can be used for collecting gold dust.

Technical Details

Dimensions (WxHxD)	250 x 150 x 150 mm (9.84 x 5.9 x 5.9 inch)
Weight	400 g

Renfert Dust Extractor Box see page 29.

Description	Art. No.
Grinding box	2625-0000