

Overview of BEGO Investment Materials

Indications and recommended liquid

Overview of BEGO investment materials

Crowns and bridges



Indication	Bellavest® SH	Bellavest® DR	Bellavest® T	BellaStar XL	Bellasan
Casting non-precious alloys	✓✓✓	✓✓✓	✓✓✓	✓✓	✓✓✓
Double crowns in non-precious alloys	✓✓✓	✓✓✓	✓✓ ¹	✓	✓✓✓
Casting precious alloys	✓✓	✓✓	✓✓	✓✓✓	✓✓
Pressable ceramics	✓✓✓	–	–	–	–
Implant prosthodontics	✓✓✓	✓✓✓	✓✓	✓✓✓	✓
CoCr partial-denture duplication with silicone	✓ ²	–	✓ ²	✓ ²	✓ ²
CoCr partial-denture hydrocolloid duplication	–	–	–	–	–

Technical data

Shock heat	✓✓✓	✓✓✓	–	✓✓✓	–
Conventional	✓✓✓	✓✓✓	✓✓✓	✓✓✓	✓✓✓
Working time ³ (20°C) [min]	4:30–5:00	5:00	5:00	3:30	7:00
Flowability [mm]	140–145	135–140	approx. 125	approx. 135	approx. 155

Scope of delivery

Accessories (contents)	REF 54257 144 × 90 g bag	REF 54861 80 × 160 g bag	REF 54202 80 × 160 g bag	REF 54362 80 × 160 g bag	REF 54270 80 × 160 g bag
	REF 70060 50 × 100 g bag				
	REF 54252 80 × 160 g bag				
	REF 54247 30 × 160 g bag				

Accessories

BegoSol® mixing liquid REF 51090 (1 liter) REF 51091 (5 liter)	–	–	✓✓✓	–	✓✓✓
BegoSol® HE mixing liquid REF 51095 (1 liter) REF 51096 (5 liter)	✓✓✓	✓✓✓	✓✓✓	–	–
BegoSol® K mixing liquid REF 51120 (1 liter) REF 51121 (5 liter)	–	–	–	✓✓✓	–
BegoSol® CC mixing liquid REF 54907 (1 liter) REF 54908 (5 liter)	–	–	–	–	–

✓✓✓ optimal · ✓✓ recommended · ✓ suitable

¹ with BegoSol® HE · ² lift-off procedure · ³ after mixing · ⁴ only conventional

Partial dentures

3D-CAD/Cast®-frames



WiroFine	Wiroplus® S	Wirovest®	Wirovest ⁺ plus	VarseoVest P plus	VarseoVest C&B
–	–	–	–	✓✓	✓✓✓
–	–	–	–	–	–
✓	✓✓	✓	✓	–	✓✓
–	–	–	–	–	–
–	–	–	–	–	✓✓
✓✓✓	✓✓✓	✓✓	✓✓	✓	–
✓✓	–	✓✓✓	✓✓✓	–	–
✓✓✓	–	–	–	✓✓✓	✓✓✓
✓✓✓	✓✓✓	✓✓✓	✓✓✓	–	✓✓
3:30	4:00	3:00	3:15	4:40	3:15
approx. 140	approx. 130	approx. 115	approx. 120	approx. 145	approx. 140

REF 54348 30 × 200g bag	REF 50248 45 × 400g bag	REF 51046 45 × 400g bag	REF 54821 45 × 400g bag	REF 54910 72 × 250g bag	REF 54894 80 × 160g bag
REF 54345 45 × 400g bag				REF 54911 60 × 300g bag	
REF 54344 15 × 400g bag				REF 54912 20 × 300g bag	

✓✓✓ ⁴	✓✓✓	✓✓✓	✓✓✓	–	–
–	–	–	–	–	–
✓✓✓	–	–	–	✓✓✓	–
–	–	–	–	–	✓✓✓

Pictures and illustrations are exemplary. Colors, symbols, design, and information on the labels and/or packaging shown may differ from reality.