

FLI® Signature Clear Bracket

FLI signature
CLEAR

RMO®'s FLI® Signature Clear Brackets offer natural aesthetic appeal to patients with the same reliable functionality preferred by doctors.



Ligation

Tie wings incorporate undercuts and downdraft for easy and secure ligation

Reduced Friction

Smooth rounded flared archslot lead-ins for reduced friction and improved sliding mechanics

Finish

Fully transparent for unmatched aesthetics without glare or reflections

Finish

Smooth rounded surface is fully polished for exceptional patient comfort

Base

Patented mushroom-style base rails deliver superior mechanical bonding with reliable and consistent debonding

Contoured base design provides optimum adaptation to tooth surfaces

Material: Polycrystalline Clear Ceramic – 99.9% Pure Alumina Oxide


Prescriptions: Roth, Bioprogressive®, MBT™



Features	Benefits
Tie wings with deep undercuts and downdraft	Easy, more secure ligation
Fully transparent	Unmatched aesthetics and true natural appearance
Smooth, fully polished surfaces	No cloudy appearance as compared to traditional ceramics
Low profile with rounded contours	Exceptional patient comfort
Smooth, flared archslot	Improved sliding mechanics – less binding
	Reduced friction over traditional ceramics
Patented mushroom-style base	Superior bond performance and strength
	Consistent debonding – greater confidence without sacrificing strength
Contoured base design	Optimum adaptation to tooth surfaces
Polycrystalline Ceramic	Improved fracture strength

FLI® Signature Clear


ROTH System Prescription



Maxillary	Torque	Angle	Offset	Hook	R / L	.018	.022
Central	12°	5°	0°	-	R	Q07200	Q07300
				-	L	Q07201	Q07301
Lateral	8°	9°	0°	-	R	Q07202	Q07302
				-	L	Q07203	Q07303
Cuspid	0°	11°	4°MO	-	R	Q07204	Q07304
				-	L	Q07205	Q07305
				HK	R	Q07206	Q07306
				HK	L	Q07207	Q07307
1st / 2nd Universal Bicuspid	-7°	0°	0°	-	R / L	Q07208	Q07308
				HK	R / L	Q07210	Q07310
Mandibular	Torque	Angle	Offset	Hook	R / L	.018	.022
Central / Lateral	0°	0°	0°	-	R / L	Q07220	Q07320
Cuspid	-11°	5°	2°MO	-	R	Q07224	Q07324
				-	L	Q07225	Q07325
				HK	R	Q07226	Q07326
				HK	L	Q07227	Q07327
1st Bicuspid	-17°	0°	0°	-	R / L	Q07228	Q07328
				HK	R / L	Q07230	Q07330
2nd Bicuspid	-22°	0°	0°	-	R / L	Q07232	Q07332
				HK	R / L	Q07234	Q07334

FLI® Signature Clear

RICKETTS® BIOPROGRESSIVE® System Prescription




Maxillary	Torque	Angle	Offset	Hook	R / L	.018
Central	22°	0°	0°	-	R / L	Q07600
Lateral	14°	8°	0°	-	R	Q07602
				-	L	Q07603
Cuspid	7°	5°	0°	-	R	Q07604
				-	L	Q07605
				HK	R	Q07606
				HK	L	Q07607
1st / 2nd Universal Bicuspid	0°	0°	0°	HK	R / L	Q07610
Mandibular	Torque	Angle	Offset	Hook	R / L	.018
Central / Lateral	0°	0°	0°	-	R / L	Q07620
Cuspid	7°	5°	0°	-	R	Q07624
				-	L	Q07625
				HK	R	Q07626
				HK	L	Q07627
1st / 2nd Universal Bicuspid	0°	0°	0°	HK	R / L	Q07610



FLI® Signature Clear

MBT™ System Prescription



Maxillary	Torque	Angle	Offset	Hook	R / L	.018	.022
Central	17°	4°	0°	-	R	Q07400	Q07500
				-	L	Q07401	Q07501
Lateral	10°	8°	0°	-	R	Q07402	Q07502
				-	L	Q07403	Q07503
Cuspid	0°	8°	0°	-	R	Q07404	Q07504
				-	L	Q07405	Q07505
				HK	R	Q07406	Q07506
				HK	L	Q07407	Q07507
	-7°	8°	0°	-	R	Q07414	Q07514
				-	L	Q07415	Q07515
				HK	R	Q07416	Q07516
				HK	L	Q07417	Q07517
1st / 2nd Universal Bicuspid	-7°	0°	0°	-	R / L	Q07408	Q07508
				HK	R / L	Q07410	Q07510
Mandibular	Torque	Angle	Offset	Hook	R / L	.018	.022
Central / Lateral	-6°	0°	0°	-	R / L	Q07422	Q07522
Cuspid	0°	3°	0°	-	R	Q07424	Q07524
				-	L	Q07425	Q07525
				HK	R	Q07426	Q07526
				HK	L	Q07427	Q07527
	-6°	3°	0°	-	R	Q07436	Q07536
				-	L	Q07437	Q07537
				HK	R	Q07438	Q07538
				HK	L	Q07439	Q07539
1st Bicuspid	-12°	0°	0°	-	R / L	Q07428	Q07528
				HK	R / L	Q07430	Q07530
2nd Bicuspid	-17°	0°	0°	-	R / L	Q07432	Q07532
				HK	R / L	Q07434	Q07534

Roth Kits	0.018	0.022
Single Case, 20 Brackets Roth FLI® Signature Clear without hooks, 5-5 U/L	KQ07200	KQ07300
Single Case, 20 Brackets Roth FLI® Signature Clear with hook on cuspid, 5-5 U/L	KQ07201	KQ07301
Single Case, 20 Brackets Roth FLI® Signature Clear without hook on cuspid / bicuspid, 5-5 U/L	KQ07202	KQ07302
Single Case, 6 Brackets Roth FLI® Signature Clear with hook on cuspid, 3-3 Upper	KQ07203	KQ07303
Single Case, 12 Brackets Roth FLI® Signature Clear without hook, 3-3 U/L	KQ07204	KQ07304
Single Case, 12 Brackets Roth FLI® Signature Clear with hook on cuspid, 3-3 U/L	KQ07205	KQ07305

Ricketts® Bioprogressive® Kits	0.018
Single Case, 20 Brackets Ricketts® Bioprogressive® FLI® Signature Clear with hook on cuspid / bicuspid, 5-5 U/L	KQ07602
Single Case, 6 Brackets Ricketts® Bioprogressive® FLI® Signature Clear with hook on cuspid, 3-3 Upper	KQ07603

MBT™ Kits	0.018	0.022
Single Case, 20 Brackets MBT™ FLI® Signature Clear without hooks, 5-5 U/L	KQ07400	KQ07500
Single Case, 20 Brackets MBT™ FLI® Signature Clear with hook on cuspid, 5-5 U/L	KQ07401	KQ07501
Single Case, 20 Brackets MBT™ FLI® Signature Clear with hook on cuspid / bicuspid, 5-5 U/L	KQ07402	KQ07502
Single Case, 6 Brackets MBT™ FLI® Signature Clear with hook on cuspid, 3-3 Upper	KQ07403	KQ07503
Single Case, 12 Brackets MBT™ FLI® Signature Clear without hook, 3-3 U/L	KQ07404	KQ07504
Single Case, 12 Brackets MBT™ FLI® Signature Clear with hook on cuspid, 3-3 U/L	KQ07405	KQ07505
Single Case, 20 Brackets MBT™ FLI® Signature Clear (Extraction Case) with 0°T cuspid, without hooks, 5-5 U/L	KQ07406	KQ07506
Single Case, 20 Brackets MBT™ FLI® Signature Clear (Extraction Case) with 0°T cuspid, hook on cuspid, 5-5 U/L	KQ07407	KQ07507
Single Case, 20 Brackets MBT™ FLI® Signature Clear (Extraction Case) with 0°T cuspid, hook on cuspid / bicuspid, 5-5 U/L	KQ07408	KQ07508
Single Case, 6 Brackets MBT™ FLI® Signature Clear (Extraction Case) with 0°T cuspid, hook on cuspid, 3-3 Upper	KQ07409	KQ07509
Single Case, 12 Brackets MBT™ FLI® Signature Clear (Extraction Case) with 0°T cuspid, without hooks, 3-3 U/L	KQ07410	KQ07510
Single Case, 12 Brackets MBT™ FLI® Signature Clear (Extraction Case) with 0°T cuspid, hook on cuspid, 3-3 U/L	KQ07411	KQ07511

