


Unlimited cure depth and no leakage propels

bulKEZ™

WITH INTELLITEK™ TECHNOLOGY


**into a class of
its own.**



**Unlimited Depth
of Cure**



**Superb
Polishability**



**Direct, Controlled
Shrinkage**



**Compatible with
a Wide Variety of
Bonding Agents**



Introducing **bulkeZ**TM with IntelliTekTM Curing Technology

The Optimal Bulk Filling Solution for Restorative Procedures

- Bulk EZTM is an easy-to-place, dual-cure composite that combines flowable cavity adaptation with high strength and wear resistance—all in one, simple step.
- Its self-cure, patent-pending IntelliTek Technology is designed to specifically control and direct shrinkage while eliminating leakage in all posterior restorations.
- Bulk EZ has reinvented bulk-fill composites by addressing the limitations in bulk fills that plague current products, such as gap formation and stress from rapid curing.


“Bulk EZ is a game changer that I believe in and endorse.”

Raymond L. Bertolotti, DDS




Unlimited Depth of Cures

- Bulk EZ ensures complete curing in the deepest cavities.
- It's designed to offer an unlimited depth of cure in one simple step without multiple layering.
- Bulk EZ's single-step process can save up to 60% of chair time over traditional light curing.




Watch the whole video by visiting
danvillematerials.com/bulkez



Well Controlled Stress and Shrinkage

- Bulk EZ™'s self-cure method offers more control over shrinkage than light cure bulk fills.
- Shrinkage is directed toward the bottom of the restoration, virtually eliminating voids and cracks.
- Bulk EZ IntelliTek™ Technology results in outstanding interproximal marginal integrity.




BULK EZ

TRADITIONAL




Excellent Physical Properties

- Bulk EZ has a flexural strength of 138MPa—among the best in dental composites.
- It has a high radiopacity of 320% AI for easy viewing on X-rays.



Compatible with a Wide Variety of Bonding Agents

- Bulk EZ's self-cure, patent-pending IntelliTek Technology ensures full compatibility with a wide variety of bonding systems.



CLASS I
PREPARED CAVITY



Superb Polishability

- In a recent study, Bulk EZ received high ratings for polishability after repeated wear cycle testing.
- Comes in 3 shades: A1, A2, and A3



CLASS I
DIRECT RESTORATION WITH BULK EZ

bulkeZ™

WITH INTELLITEK™ TECHNOLOGY

Experience a bulk-fill composite in a class of its own. Call **800.827.7940** or visit **Danvillematerials.com/bulkez**

Bulk EZ™ Ordering Information:

94778 Bulk EZ Trial Kit

Bulk EZ A1 syringe (x1), Bulk EZ A2 syringe (x1), Bulk EZ A3 syringe (x1), Bulk EZ 17 gauge mixing tips (x10), Bulk EZ 19 gauge mixing tips (x10).

94775 Bulk EZ Refill A1

Bulk EZ A1 syringe (x2), Bulk EZ 17 gauge mixing tips (x6), Bulk EZ 19 gauge mixing tips (x6).

94776 Bulk EZ Refill A2

Bulk EZ A2 syringe (x2), Bulk EZ 17 gauge mixing tips (x6), Bulk EZ 19 gauge mixing tips (x6).

94777 Bulk EZ Refill A3

Bulk EZ A3 syringe (x2), Bulk EZ 17 gauge mixing tips (x6), Bulk EZ 19 gauge mixing tips (x6).

94792 Bulk EZ Refill 17

Gauge mixing tips (x25)

94793 Bulk EZ Refill 19

Gauge mixing tips (x25)

Bulk EZ™ Accessory Products:



Mega V™ Contact Matrix Ring

Danville's new Mega V Matrix system is the perfect solution for your Class II restorations. The Mega V Ring, with its unique inclined hoop and superior strength, allows you to perform Class II restorations both quickly and efficiently.



Prelude One™ Universal Adhesive

Prelude One is a light curable single component adhesive used for bonding to all dental surfaces. Its extremely high bond strengths and compatibility with dual-cure composites (without requiring an additional activator) make it an ideal bonding agent.

ZD ZEST DENTAL SOLUTIONS®

ZEST | DANVILLE MATERIALS | PERIOSCOPY

www.danvillematerials.com | Phone: 800.827.7940 | Fax: 925.973.0764