



FenderMate® D-Coat

Prod No: 602790	FenderMate D-Coat, refill, dark green, 18 pcs
Prod No: 602791	FenderMate D-Coat, refill, light green, 18 pcs
Prod No: 602792	FenderMate D-Coat, refill, dark blue, 18 pcs
Prod No: 602793	FenderMate D-Coat, refill, light blue, 18 pcs
Prod No: 602794	FenderMate D-Coat, assortment, 4 x 18 pcs

FenderMate D-Coat is a new addition to the FenderMate series. Our unique no-ring matrix system is now redesigned with a multilayer diamond coating combining a better visual contrast with anti-stick properties. Your high-quality class II composite restorations can now be achieved without adhesion to the matrix. The dark matt color will facilitate your work by giving excellent contrast. Like the original FenderMate, the pre-curved matrix includes a contact point and comes attached to a specially designed wedge which allows tight cervical sealings with no overhang.

FenderMate D-Coat is available in refill and assortments. The four colors allow an easy distinction between left or right curve (green/blue), and wedge size (light = thin, dark = thick).

4

Reasons to use FenderMate® D-Coat

Non-sticking coating

The amorphous diamond coating prevents adhesives from sticking to the surface. Pulling out the matrix after light curing will be an effortless task.

Dark & matt surface

The dark and matt coating contrasts with the color of the composite resin and reduces light reflection. Clinicians can therefore perform restorations under excellent visual conditions.

One piece design

Matrix and wedge combined for easier and safer placement.

No cervical overhang

The wing design helps maintain proximal separation and secures superior cervical adaptations.
