

**Zekrya**

**E 0151**

**Applications:**

- Sectioning an impacted wisdom tooth prior to its extraction
- Separating roots
- Removing a broken root stump
- Cutting a tooth horizontally at neck level
- Reworking temporary teeth of acrylic material
- Preparing a shoulder on the vestibular side of the tooth

**Caractéristiques:**

- Active laterally (no vibrations on account of its 6 helical cutting edges)
- Active at end by means of its hemispherical, cutting tip



6 helical cutting edges: active laterally

Hemispherical, cutting tip

