

Imprint™ 4 Bite

VPS Bite Registration Material



As fast and accurate
as it can be

Imprint™ 4 Bite VPS Material offers you:

- ▶ **Good flow behaviour**
For high accuracy
- ▶ **Fast intra-oral setting**
For a faster procedure
- ▶ **High final hardness and low flexibility**
For easy cutting and trimming and accurate stone model positioning
- ▶ **Fresh material colour and scannability**
For good readability and easy access to the digital workflow
- ▶ **Cartridge delivery**
For simple and convenient mixing and dispensing

Developed for highly accurate bite registrations.

With its great flow properties and a high enough viscosity, Imprint 4 Bite material does not flow deeply into interdental areas and delivers highly accurate bite registrations. A high final hardness and low flexibility after setting allow for easy cutting and trimming as well as accurate positioning of stone models.

Saves time with a fast intra-oral setting.

Imprint 4 Bite offers a pleasant working time and a short setting time of only 60 seconds, keeping procedures moving quickly while enhancing patient comfort. And, it helps to reduce inaccuracies due to jaw movements during setting.

	Mixing and application	Setting time from start of mixing*
Imprint™ 4 Bite VPS Bite Registration Material	0:20	1:00

*The stated times are valid at 23°C/74°F. Times depend on ambient temperature.

Indications:

Universal VPS bite registration material – ideally suited for occlusal recording including optical registration of occlusal data for CAD/CAM systems.



Scan of bite registration made with Imprint™ 4 Bite VPS Material.

Get full control of your impression work.

3M ESPE impression solutions are designed to add speed, convenience and precision to the entire procedure.



Imprint™ 4 Bite VPS Material is part of the successful family of Imprint 4 VPS impression materials. **Try also the versatile range of Imprint™ 4 VPS materials for your everyday precision impressions as well as Imprint 4 Preliminary VPS Material.**

Ordering Information



Imprint™ 4 Bite Refill (71529)

Item No.	Product Information
71529	Imprint™ 4 Bite Refill 2 cartridges – 50ml each, 5 Garant™ Mixing Tips green, 5 Garant™ Intra-oral Tips Bite
71460	Garant intraoral tip Bite

www.3MESPE.com

3M ESPE

3M Deutschland GmbH
Location Seefeld
3M ESPE · ESPE Platz
82229 Seefeld · Germany
info3mespe@mmm.com
www.3MESPE.com

3M, ESPE, Garant and Imprint are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada.

© 2014, 3M. All rights reserved.

XX (0X.2014)