

PLANMECA DENTAL UNITS

HYGIENIC COVER FOR PANASONIC INTRAORAL CAMERA PART NUMBER 10037205

INSTRUCTIONS FOR USE

The hygienic cover has been manufactured from 2 plastic films.

The plastic film closest to the paper sheet is the optically good one.

Always place the camera so that the lens of the camera points toward the paper sheet (as demonstrated in the picture below).



Remove the paper sheet before using the intraoral camera.

The order code for the hygienic cover is 10037205.

EN

Camera Sleeve



A protective sleeve for intra-oral cameras.
The product may only be used as a protective sleeve for intra-oral cameras.
The manufacturer or its resellers cannot be held liable for any other use.
Store at room temperature, protected from dust and humidity. Do not expose to direct sunlight. Reuse will lead to contamination and possible spread of infection.

Camera Sleeve 305, 500 pcs

CE REF 413190305